HUNTLEIGH

D900/D920/D930

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Mode d'emploi

Bruksanvisning

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nwijzing

INSTRUCTIONS FOR USE

Bruksan

matları

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





WARNING/CAUTIONS & SAFETY

Do not use Dopplers in the presence of flammable gases such as anaesthetic agents.

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Do Not :

- Immerse in any liquid, (except D920/D930 probe)
- use solvent cleaner,
- use high temperature sterilising processes (such as autoclaving),
- use E-beam or gamma radiation sterilisation.

This product contains sensitive electronics, therefore, strong radio frequency fields could possibly interfere with it. This will be indicated by unusual sounds from the loudspeaker. We recommend that the source of interference is identified and eliminated.

Do not use on the eye.

If outputs 1 and 2 are used simultaneously, the system should comply with EN60601-1-1. Any equipment connected to outputs 1 or 2 should comply with EN60601-1, EN60950, EN60065, EN60335 or EN61010.

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Do not attempt to recharge normal dry-cell batteries. They may leak, cause a fire or even explode. Dopplex Dopplers are screening tools to aid the healthcare professional and should not be used in place of normal vascular or fetal monitoring. If there is doubt as to vascularity or fetal well-being after using the unit, further investigations should be undertaken immediately using alternative techniques. We recommend that exposure to ultrasound should be kept As Low As Reasonably Achievable - (ALARA guidelines). This is considered to be good practice and should be observed at all times.

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Means Attention consult this manual. Refer to safety section.

GENERAL OPERATION & MAINTENANCE

Refer to diagram on inside front cover for Doppler Measuring sites and Recommended Probes.

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To replace battery see diagram.

After Use

1. Press and release the On/Off button. If you forget to switch the unit off, it will automatically shut-off after 5 minutes.
2. Refer to the cleaning section before storing or using the unit on another patient.
3. Store unit together with probe and accessories in the soft carry case provided.

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Cleaning

Remove excess gel before parking the probe. Clean with a damp cloth impregnated with mild detergent. Do not allow fluid to seep into either unit.

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Coupling Gel

Use water based ultrasound gel ONLY.

OPERATING INSTRUCTIONS VASCULAR

Clinical Use


Apply a liberal amount of gel on the site to be examined. Place the probe at 45° to the skin surface over the vessel to be examined. Adjust the position of the probe to obtain the loudest audio signal. High pitched pulsatile sounds are emitted from arteries while veins emit a non-pulsatile sound similar to a rushing wind. For best results, keep the probe as still as possible once the optimum position has been found. Adjust the audio volume as required.

Waveform Recording (Mini only)

The Mini Dopplex can be connected to a chart recorder by using a suitable lead that plugs into the Waveform socket. (For more information on the chart recorder, refer to the manufacturer's handbook).

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Two levels of calibration pulses are provided. While Cal is active, calibration tones synchronised to the calibration pulses are superimposed on the audio signal.

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VP4HS	4MHz $\pm$ 1% for deep lying vessels
VP5HS	5MHz $\pm$ 1% for deep lying vessels and oedematous limbs
VP8HS	8MHz $\pm$ 1% for peripheral vessels
VP10HS	10MHz $\pm$ 1% for specialist superficial applications.
EZ8	8MHz $\pm$ 1% "Widebeam" for peripheral vessels.

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In early pregnancy a full bladder may improve sound detection. In later pregnancy the best signals are generally located higher on the abdomen. The fetal heart sounds like a galloping horse at approximately twice the maternal rate. A wind-like sound is heard from the placenta.

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Note that early detection (9-12 weeks) is highly dependent on operator skill and other factors and may not always be possible.

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Audio Dopplex (D920/D930)

The Audio Dopplex has a 2MHz/3MHz waterproof probe for underwater use.

Probes

OP2HS	2MHz $\pm$ 1%
OP3HS	3MHz $\pm$ 1%

WARRANTY & SERVICE

Huntleigh Healthcare's standard terms and conditions apply to all sales. A copy is available on request.

These contain full details of warranty terms and do not limit the statutory rights of the consumer.

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1. Clean the product, as described in the cleaning section.
2. Pack it in suitable packing.
3. Attach the decontamination certificate (or other written statement declaring that the product has been cleaned) to the outside of the package.
4. Mark the package "Service Department - Mini/Audio Dopplex".

UK

For further details refer to the NHS document HSG(93) 26.

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Or your local distributor.

Manufactured in the UK by Huntleigh Healthcare Ltd. As part of the ongoing development programme the company reserves the right to modify specifications and materials of the Mini Dopplex or Audio Dopplex without notice.

Dopplex, Huntleigh and 'H' logo are registered trademarks of Huntleigh Technology Ltd 2004.

© Huntleigh Healthcare Ltd. 2004

TECHNICAL SPECIFICATION

Product Name:	Mini/Audio Dopplex®		
Model No.:	D900/D920/D930		
Physical Characteristics:	Weight : 295 gms Height : 140mm	Width : 74mm	Depth : 27mm
Max. Audio Output:	500mW rms typical		
Max Headphones Output:	25 mW rms (32Ω headphones) (max. applied voltage +9Vdc)		
Auto Shut-Off:	After 5 minutes continuous operation		
Waveform Outputs: (Mini Dopplex only)	Zero crosser, 0.5V/kHz, 3.5V full scale Cal levels:1kHz + 2kHz, (max. applied voltage +9Vdc) 3.5mm mono socket		
Battery Type Recommended:	9 volt alkaline manganese (e.g. MN1604)		
Battery Life:	Typically, 500 x 1 minute examinations		
Audio Dopplex Probe:	Waterproof to IPX7		
Complies With:	BS5724: Part 1 : 1989, IEC 601-1: 1988, EN60601-1		
IEC601-1 Classification:	Type of shock protection	Internally powered equipment	
	Degree of shock protection	Type B equipment	
	Protection against water ingress	Ordinary equipment except D920/D930 probe: IPX7 - protected against immersion	
	Degree of safety in presence of flammable gases	Equipment not suitable for use in presence of flammable gases	
	Mode of operation	Continuous	
Operating Temperature:	+10°C to +30°C		
Storage Temperature:	-10°C to +40°C		

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ENGLISH
(USA)

WARNING!

Federal law restricts this device to sale by or on the order of a licensed practitioner.

IMPORTANT

Before using your Audio/Mini/Flexi Dopplex, please study this manual carefully and familiarize yourself with the controls, displays, features and operating techniques.

Do not use Dopplers in the presence of flammable gases such as anesthetic agents.

Do not use in the sterile field unless additional barrier precautions are taken.

Do not :

- Immerse in any liquid (Except **D920/D930** probe)
- Use solvent cleaner
- Use high temperature sterilizing processes (such as autoclaving)
- Use E-beam or gamma radiation sterilization

This product contains sensitive electronics, therefore, strong radio frequency fields could possibly interfere with it. This will be indicated by unusual sounds from the loudspeaker. We recommend that the source of interference be identified and eliminated.

Do not use on the eyes or scrotum.

If outputs 1 and 2 are used simultaneously, the system should comply with EN60601-1-1.

Any equipment connected to outputs 1 or 2 should comply with EN60601-1, EN60950, EN60065, EN60335 or EN61010.

Do not dispose of batteries in fire as this can cause them to explode.

Do not attempt to recharge normal dry-cell batteries. They may leak, cause a fire or even explode.

Dopplex Dopplers are screening tools to aid the healthcare professional and should not be used in place of normal vascular or fetal monitoring. If there is doubt as to vascularity or fetal well-being after using the unit, further investigations should be undertaken immediately using alternative techniques.

Audio Dopplex (D920/D930)

The main unit is not waterproof and must not be immersed. For underwater use where contamination or cross-infection may occur, additional barrier precautions must be taken.

Acoustic Safety

Continuous wave Doppler ultrasound instruments such as the D900, D920 & D930 have been used extensively for medical diagnosis in the United States for over 25 years. Throughout this period, there have been no reports of adverse effects to patients or instrument operators at the acoustic intensities recommended for diagnostic use. Despite this highly favorable safety experience, available data are not conclusive and the possibility remains that unwanted biological effects might be identified in the future.

Authorities therefore recommend that ultrasound procedures be performed in accordance with the "ALARA" principle, which states that the energy delivered to the patient should always be kept As Low As Reasonably Achievable. With the D900, D920 & D930, the transmitted acoustic power is fixed and cannot be adjusted by the operator. Therefore, the user can best observe the ALARA principle by ensuring that each examination is medically indicated and by limiting the duration of the study to the extent appropriate for the clinical objectives.

WARNING/CAUTIONS & SAFETY

Acoustic Safety (continued)

Acoustic intensity data ($I_{SPTA,3}$) for probes available for use with the D900, D920 & D930 are summarized in the following table. The values cited are based on measurements in water using a calibrated hydrophone and are stated as the estimated derated intensities. The derated intensity constitutes the most biologically relevant parameter available since true determinations of actual absorbed dose in tissue would require invasive measurement techniques. The derated intensity is therefore calculated mathematically using a derating factor consisting of a constant (the assumed attenuation coefficient) and allowing for the frequency of the probe and the distance from the probe face to the hydrophone. The calculated derated intensity values for the D900, D920 & D930 compare very favorably with previously reported acoustic safety data for Doppler ultrasound instruments and are appropriate for all clinical applications recommended in this manual. As the operating mode of the Doppler range of probes is continuous wave, I_{SPPA} figures are not applicable.

Acoustic Output Table, Track1, Non-Auto-Scanning Mode						
Model	Max. Value $I_{SPTA,3}$	W_o	f_c	Z_{sp}	A-6, (Z_{sp})	EBD
OP2 HS	55	41	2.0	2.5	1.1	1.2 x 2.5
OP3 HS	55	32	3.0	2.0	1.2	1.2 x 2.5
VP4 HS	92	7.5	4.0	0.8	0.14	0.365 x 0.8
VP5 HS	92	8.2	5.0	0.8	0.12	0.365 x 0.8
VP8 HS	92	4.0	8.0	0.48	0.026	0.215 x 0.5
EZ8	92	14.3	8.0	0.67	0.064	0.635 x 0.22
VP10 HS	92	1.4	10.0	0.48	0.022	0.215 x 0.5

NOTES

- Measurement uncertainty: varies with probe and measurement
Random - typically $\pm 20\%$ (max. $\pm 32\%$)
Systematic - typically $\pm 6.5\%$ (max. $\pm 8\%$)

Definition of Terms

- $I_{SPTA,3}$** is the derated spatial-peak, temporal-average intensity (milliwatts per square centimeter)
- W_o** is the ultrasonic power (milliwatts)
- f_c** is the center frequency (MegaHertz)
- Z_{sp}** is the axial distance used to calculate the derated intensity (centimeters)
- A-6 (Z_{sp})** is $(\pi/4) \times (X-6 \times Y-6)$ where X-6, Y-6 are respectively the in-plane (azimuthal) and out-of-plane (elevational) -6dB dimensions in the X-Y plane where Z_{sp} is found (centimeters)
- EBD** are the entrance beam dimensions for the azimuthal and elevational planes (centimeters)

GENERAL OPERATION & MAINTENANCE

ENGLISH
(USA)

Refer to diagram on Page 2 for Doppler Measuring Sites and Recommended Probes.

Indicators

Green LED - indicates power ON.

Yellow LED - flashes when battery is low.

During Use

An automatic noise reduction feature operates on low level signals to improve sound quality.

To replace battery, see diagram on page 2.

After Use

1. Press and release the On/Off button. If you forget to switch the unit off, it will automatically switch off after 5 minutes.
2. Refer to the cleaning section before storing or using the unit on another patient.
3. Store unit together with probe and accessories in the soft carry case provided.

The control unit and the body of the probe are robust and require no special handling. However, the probe tip is delicate and must be handled with care.

Cleaning

Remove excess gel before replacing the probe. Clean with a damp cloth impregnated with mild detergent. Do not allow fluid to seep into either unit.

For disinfection use a soft cloth with Sodium Hypochlorite 1000ppm or alcohol.

Alcohol saturated swabs can be used on the probe tip.

Please be sure to check your local control of infection policies, or any equipment cleaning procedures.

Phenolic, detergent based disinfectants containing cationic surfactants, ammonia based compounds, or antiseptic solutions such as Steriscol or Hibiscrub should never be used on any part of the system.

Coupling Gel

Use water-based ultrasound gel ONLY.

D920/D930

Gel not required when probe used under water.

**Attention! Consult User Manual.
Refer to Safety Section.**

OPERATING INSTRUCTIONS VASCULAR

Clinical Use

Apply a liberal amount of gel on the site to be examined. Place the probe at 45° to the skin surface over the vessel to be examined.

Adjust the position of the probe to obtain the loudest audio signal. High pitched pulsatile sounds are emitted from arteries while veins emit a non-pulsatile sound similar to a rushing wind.


For best results, keep the probe as still as possible once the optimum position has been found. Adjust the audio volume as required.

Waveform Recording (D900 Only)

The D900 can be connected to a chart recorder by using a suitable lead that plugs into the Waveform socket. (For more information on the chart recorder, refer to the manufacturer's handbook).

Cal Function

The baseline and sensitivity of the chart recorder can be set up using the Cal function. This generates a zero velocity baseline and a sequence of pulses as shown below:

Two levels of calibration pulses are provided. While Cal is active, calibration tones synchronized to the calibration pulses are superimposed on the audio signal.

If the button is pressed and held, then the calibration sequence is repeated. When the calibration sequence is complete, normal operation is resumed.

OPERATING INSTRUCTIONS VASCULAR

Probes

Five probes are available for vascular examinations:

VP4 HS	4MHz $\pm$ 1% for deep lying vessels
VP5 HS	5MHz $\pm$ 1% for deep lying vessels and edematous limbs
VP8 HS	8MHz $\pm$ 1% for peripheral vessels
VP10 HS	10MHz $\pm$ 1% for specialist superficial applications.
EZ8	8MHz $\pm$ 1% "Widebeam" for peripheral vessels.

OPERATING INSTRUCTIONS OBSTETRICS

Clinical Use

Apply a liberal amount of gel to the abdomen. Place the faceplate of the probe flat against the abdomen above the symphysis pubis. Adjust the probe to obtain an optimum audio signal ideally by angling the probe around. Avoid sliding it over the skin.

In early pregnancy, a full bladder may improve sound detection. In later pregnancy the best signals are generally located higher on the abdomen. The fetal heart sounds like a galloping horse at approximately twice the maternal rate. A wind-like sound is heard from the placenta.

The Dopplex units can detect fetal heart sounds from as early as 9 to 10 weeks gestation.

Note that early detection (9 – 12 weeks) is highly dependent on operator skill and other factors and may not always be possible.

Mini Dopplex only

Obstetric mode is automatically selected when an obstetric probe (**OP2HS** or **OP3HS**) is connected to the control unit.

In this mode the waveform output and Cal function are disabled.

Audio/Flexi Dopplex (D920/D930)

Both the Audio Dopplex and the Flexi Dopplex have waterproof probes for underwater use.

Probes for obstetric applications

OP2 HS : 2MHz $\pm$ 1%
OP3 HS : 3MHz $\pm$ 1%

WARRANTY

a) HUNTLEIGH HEALTHCARE LLC. HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES (INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE) AND ANY AGREEMENTS, REPRESENTATIONS, AFFIRMATIONS, OR WARRANTIES, WHETHER ORAL OR WRITTEN, MADE BY ANY AGENT, EMPLOYEE OR REPRESENTATIVE OF HUNTLEIGH HEALTHCARE INC., UNLESS SPECIFICALLY SET FORTH IN THIS PARAGRAPH. HUNTLEIGH HEALTHCARE LLC. SHALL NOT BE LIABLE FOR BREACH OF CONTRACT ARISING FROM ANY DEFECT IN MATERIAL OR WORKMANSHIP OF THE GOODS. ALL LEGISLATION RELATING TO EXPRESS AND IMPLIED WARRANTIES OR OTHER OBLIGATIONS ON THE PART OF HUNTLEIGH HEALTHCARE LLC. THAT MAY BE LAWFULLY EXCLUDED ARE HEREBY EXCLUDED.

b) Notwithstanding the foregoing, Huntleigh Healthcare Llc.'s sole warranty is that the Goods shall be free from defects in material and workmanship for a period of three (3) years (excluding probe head and retractile cable which are warranted for one (1) year), following delivery of such Goods to the original purchaser; **provided** that the Goods were used in an appropriate and reasonable manner during such period and **provided further** that Huntleigh Healthcare Llc. shall in no event be liable to Customer for defective Goods if: (i) the Goods are damaged in the course of shipping; (ii) any defect is caused wholly or to any material extent by Customer's negligence, misuse, failure to use the Goods properly or use of the Goods in conjunction with any accessory not approved for use with the Goods by Huntleigh Healthcare Llc.; (iii) the Goods are damaged as a result of improper maintenance, failure to follow manufacturer's instructions, including without limitation those on washing and cleaning, or failure to follow necessary routine maintenance procedures; or (iv) the Goods are altered, repaired or dismantled other than with manufacturer's written authorization using its approved procedures or by any party other than manufacturer's properly qualified and trained technicians.

c) Customer must provide written notice to Huntleigh Healthcare Llc. within said warranty period of any defect in the Goods. Upon Huntleigh Healthcare Llc.'s written request, Customer must return such Goods adequately packed (in their original packing) and fully insured to Huntleigh Healthcare Llc.'s place of business and shall be responsible for all shipping costs incurred therein.

Customer's exclusive remedy and Huntleigh Healthcare Llc.'s exclusive liability for any claim for loss, damage or destruction resulting from any defects in materials and workmanship shall be limited to repair, service, adjustment or replacement (at Huntleigh Healthcare Llc.'s option) of any nonconforming or defective Goods. Huntleigh Healthcare Llc. will have a reasonable time to repair, service or replace such Goods. Any Goods returned to Huntleigh Healthcare Llc. which are found not to be defective in breach of the warranty in Subsection (b) above, shall be returned to Customer in the manner described in this subsection.

d) IN NO EVENT SHALL HUNTLEIGH HEALTHCARE LLC. BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSSES OR DAMAGES (INCLUDING BUT NOT LIMITED TO ECONOMIC LOSS, LOSS OF PROFITS OR SPECIAL DAMAGES) ARISING OUT OF OR INCURRED BY CUSTOMER IN CONNECTION WITH THE PURCHASE OF HUNTLEIGH HEALTHCARE LLC.'S GOODS EVEN IF HUNTLEIGH HEALTHCARE LLC. HAS BEEN ADVISED OR HAS KNOWLEDGE OF THE POSSIBILITY OR EXTENT OF SUCH DAMAGES SUFFERED OR INCURRED BY CUSTOMER OR ANY END USER AS A RESULT OF OR IN CONNECTION WITH ANY BREACH OF THESE TERMS AND CONDITIONS BY HUNTLEIGH HEALTHCARE LLC. OR ANY TORT (INCLUDING BUT NOT LIMITED TO STRICT LIABILITY OR NEGLIGENCE) COMMITTED BY HUNTLEIGH HEALTHCARE LLC., ITS AGENTS OR REPRESENTATIVES IN CONNECTION WITH THESE TERMS AND CONDITIONS OR ANY CONTRACT WITH CUSTOMER FOR THE SUPPLY OF GOODS.

WARRANTY & SERVICE

e) Customer shall not create, directly or indirectly, any warranty obligations on the part of Huntleigh Healthcare Inc. to the customers of Customer, and in particular, without limiting the foregoing, Customer agrees not to pass on to its customers any warranties beyond or in addition to those given by Huntleigh Healthcare Llc. to Customer hereunder. Where the Customer is a dealer in the Goods, it shall be responsible for the labor cost of all repairs and Huntleigh Healthcare Llc. shall be responsible for providing all repair parts during said three (3) year warranty period (excluding probe head and retractile cable which are warranted for one (1) year). The dealer shall provide written verification of warranty repairs including the original invoice number, date of purchase, description of repairs, name of its customer and date of sale to such customer.

f) Customer shall be deemed to have full knowledge of the nature and properties of the Goods ordered and of any hazards they involve and the proper treatment, storage and handling thereof. Any technical advice furnished by Huntleigh Healthcare Llc. or its representatives or agents is given only on the basis that it is followed at the Customer's own risk.

SERVICE

Service Returns

There are NO USER SERVICEABLE PARTS inside the control unit or probe.

If for any reason your Dopplex Unit is being returned, please:

1. Clean the product, as described in the cleaning section.
2. Pack in suitable packing.
3. Attach the decontamination certificate, (or other written statement declaring that the product has been cleaned), to the outside of the package.
4. Call Huntleigh Service Dept. in Eatontown, NJ for a return authorization number.

For service, maintenance and any questions regarding this, or any other Huntleigh Healthcare Dopplex product, please contact:

Service Department
Huntleigh Healthcare LLC
40 Christopher Way
Eatontown, NJ 07724-3327
Tel: (800) 223-1218
Local: (732) 578-9898
Fax: (732) 578-9889

Or your local distributor.

Manufactured in the UK by Huntleigh Healthcare.

As part of the ongoing development programme, the company reserves the right to modify specifications and materials of the Audio Dopplex, Mini Dopplex and Flexi Dopplex without notice.


Dopplex, Huntleigh and 'H' logo are registered trademarks of Huntleigh Technology LTD.

©Huntleigh Healthcare 2005

TECHNICAL SPECIFICATIONS

Product Name:	Mini/Audio/Flexi Dopplex®		
Model No.:	D900/D920/D930		
Physical Characteristics:	Weight :	10oz including probe and battery	
	Height :	5.5"	
	Width :	2.9"	
	Depth :	1.1"	
Max. Audio Output:	500mW rms typical		
Max Headphones Output Power:	25 mW rms max. (32Ω headphones) (max.applied voltage +5Vdc)		
Auto Shut-Off:	After 5 minutes of continuous operation.		
Waveform Outputs :	Zero crosser, 0.5V/kHz, 3.5V full scale. Cal levels: 1kHz + 2kHz, (Max. applied voltage +9Vdc) 3.5mm mono socket.		
(Mini Dopplex Only)			
Battery Type Recommended:	9 volt alkaline manganese (e.g. MN1604)		
Battery Life:	Typically, 500 x 1 minute examinations		
Audio/Flexi Dopplex Probe :	Waterproof to IPX7		
Complies With:	BS5724: Part 1 : 1989, IEC 601-1: 1988, EN60601-1		
IEC601-1 Classification:			
Type of shock protection	Internally powered equipment		
Degree of shock protection	Type B equipment		
Protection against water ingress	Ordinary equipment except D920/D930 probe: IPX7 Protected against immersion		
Degree of safety in presence of flammable gases	Equipment not suitable for use in presence of flammable gases		
Mode of operation	Continuous		
Operating Temperature:	+10°C to +30°C (50°F to 86°F)		
Storage Temperature:	-10°C to +40°C (14°F to 104°F)		

D900/D920/D930

	D900	D920/D930	ENGLISH	FRANÇAIS	DEUTSCH	ITALIANO	SVENSKA	NEDERLANDS	ESPAÑOL	DANSK	SUOMI
1	●	●	Head-phone Socket	Prise casque	Buchse Kopfhörer	Presa cuffia	Uttag för hörlurar	Hoofdstel-contact	Entrada de auriculares	Hoved-telefon holder	Kuulo-kevastake
2	●		Waveform Socket	Prise générateur d'ondes	Buchse Meßkurven	Presa segnale	Vågform-suttag	Golvform-contact	Salida para gráfica	wave form holder	Aaltomuot-ovastake
3	●	●	LED Panel	Panneau de LED	LED-Feld	Pannello del LED	LED -indikator	LED-Paneel	Panel LED	LED Panel	Valodio-dinäyttö (LED)
4	●	●	Loud-speaker	Haut-parler	Laut-sprecher	Alto-parlante	Högtalare	Luid-sprecker	Altavoz	højtaler	Kaiutin
5	●	●	On/Off Button	Bouton Marche/Arrêt	Ein/Ausschalter	Interruttore on/off	På/Av knapp	Aan/uitschakelaar	Botón de Encendido/Apagado	tænd og sluk knap	On/Off-painike
6	●	●	Volume Control	Contrôle du volume	Volume Control	Regolatore di volume	Volym-kontroll	Volume-knop	Control du Volumen	Volume kontrol	Äänenvoimakk-uussäädin
7	●		Cal Button	Bouton calibrage	CAL-Schalter	Pulsante Cal	Kalibreringsknapp	Cal-knop	Botón de Calibrado	Kalibrerings knap	Kal-painike
8	●	●	Battery Compartment	Logement pile	Batteriefach	Vano batteria	Batterifack	Batterij-holte	Compartmento para pila	Batteri rum	Paristolokero
9	●	●	Pocket Clip Attachment	Clip	Zub. Taschenclip	Fermaglio da taschino	Fästånord-ning för ficka	Zakclip	Clip para bolsillo	Lomme clips	Taskupidikeliitin
10	●	●	Probe Holder	Rangement sonde	Sondenhalterung	Reggi-sonda	Probe-parkering	Positie van sonde	Compartmento para sonda	sensor holder	Koettimen pidin

Visit our Website for
downloadable clinical &
educational material and
product information

HUNTLEIGH ...with people in mind

Huntleigh Healthcare Ltd. - Diagnostic Products Division
35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom
T: +44 (0)29 20485885 F: +44 (0)29 20492520
E: sales@huntleigh-diagnostics.co.uk
W: www.huntleigh-diagnostics.com

DISTRIBUTED IN GERMANY BY:

HNE Huntleigh Nesbit Evans Healthcare GmbH
Industriering Ost 66, 47906 Kempen
T: +49 (0) 2152-5511-10 F: +49 (0) 2152-5511-20
E: service@hne-healthcare.de

DISTRIBUTED IN THE USA BY:

Huntleigh Healthcare Inc
40 Christopher Way, Eatontown New Jersey 07724-3327
T: (800) 223-1218 F: (732) 578-9889

Registered No: 942245 England. Registered Office: 310-312
Dallow Road, Luton, Beds, LU1 1TD

© Huntleigh Healthcare Limited 2004

MEMBER OF THE GETINGE GROUP

® and ™ are trademarks of Huntleigh Technology Limited
As our policy is one of continuous improvement, we reserve the right to modify designs without
prior notice.