Fingertip Pulse Oximeter

English

General Description

Oxygen binds to hemoglobin in red blood cells when moving through the lungs. It is transported throughout the body as arterial blood. A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO₂. A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO₂ level.

Measurement Principle

Principle of the oximeter is as follows: The pulse oximeter works by applying a sensor to a pulsating arterioles vascular bed. The sensor contains a dual light source and photo detector. The one wavelength of light source is 660nm, which is red light; the other is 905nm, which is infrared-red light. Skin, bone, tissue and venous vessels normally absorb a constant amount of light over time. The photo detector in finger sensor collects and converts the light into an electronic signal which is proportional to the light intensity. The arterioles bed normally pulsates and absorbs variable amounts of light during systole and diastole, as blood volume increases and decreases. The ratio of light absorbed at systole and diastole is translated into an oxygen saturation measurement. This measurement is referred to as SpO₂.

Diagram of Operation Principle

Red and Infrared-ray Emission Tube

2. Red and Infrared-ray Reception Tube

Precautions For Use

- Before use, care must read the manual.
- Operating the fingertip pulse oximeter may be affected by the use of an electro-surgical unit (ESU).
- The fingertip pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO₂ measurement.
- Do not use the fingertip pulse oximeter in an MRI or CT environment.
- Do not use the fingertip pulse oximeter in situations where alarms are required. The device has no alarms. It is not for continuous monitoring.
- Do not use the fingertip pulse oximeter in an explosive atmosphere.
- The fingertip pulse oximeter intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
- In order to ensure correct sensor alignment and skin integrity, the maximum application time at a single site for our device should be less than half an hour.
- Do not sterilize the device using autoclaving, ethylene oxide sterilizing, or immersing the device in liquid. The device is not intended for sterilization.
- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- This device is not compatible with IEC 60601-1-2:2010 and its amendments concerning medical electrical equipment and/or systems. However, because of the provision of radio-frequency communications equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device.
- Portable and mobile RF communications equipment can affect medical electrical equipment.
- This equipment is not intended for use during patient transport outside the healthcare facility.
- This equipment should not be used adjacent to or stacked with other equipment.
- It may be unsafe to use this device in combination with other equipment not described in the instructions for use
- Interconnect the equipment with other equipment not described in the instructions for use
- Disassemble, repair or modify the equipment
- These materials that contact with the patient's skin contain medical silicone and ABS plastic enclosure are all pass the ISO10993-5 Tests for invitro cytotoxicity and ISO10993-10 Tests for irritation and delayed-type hypersensitivity.
- When the signal is not stable, the reading may inaccurate. Please do not reference.

Rx only: "Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed practitioner."

Contraindication

It is not for continuous monitoring.

Inaccurate measurements may be caused by

- Significant levels of dysfunctional hemoglobin (such as carboxy-, hemoglobin or methemoglobin).
- Intravascular dyes such as indocyanine green or methylene blue.
- High ambient light. Shield the sensor area if necessary.
- Excessive patient movement.
- High-frequency electro-surgical interference and defibrillators.
- Venous pulsations.
- Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
- The patient is in cardiac arrest or is in shock.
- Fingernail polish or false fingernails.
- Weak pulse quality (low perfusion).
- Low hemoglobin.

Product Features


- Simple to operate and convenient to carry.
- Small volume, light weight and low power consumption.
- Dual color OLED displays SpO₂, PR, PI (Perfusion Index), Pulse bar, and waveform.
- 4.7 display modes.
- Level 1-10 adjustable brightness.
- 2.5cs AAA alkaline batteries, real-time battery status indication.
- Weak or unstable signal prompt provides more accurate measurements.
- The device automatically shuts off after no operation in 8 seconds when "finger out" displays.
- Multi-patient reusability.

Intended Use

The Fingertip Pulse Oximeter is a handheld non-invasive device intended for spot-checking of oxygen saturation of arterial hemoglobin (SpO₂) and Pulse Rate of adult, adolescent and child patients in hospitals, hospital-type facilities and homecare.


Operation Instructions

- Install two AAA batteries according to the Battery Installation Instructions.
- Place the finger firmly into the rubber opening of the Pulse Oximeter.
- Press the switch button once on front panel to turn the pulse oximeter on.
- Keep your hands still for the reading. Do not shake your finger during the test. It is recommended that you do not move your body while taking a reading.
- Read the data from the display screen.
- The display modes are as follows:

Notes:
1. Short press the power button to switch the display modes.
2. Long press the power button to adjust the brightness of the oximeter. There are 10 levels of brightness. The default is level four.
3. Take out your finger, the screen displays "FingerOut". It means the measurement ending.

Reading your results

Note:
1. The pulse bar less than 30% indicates signal inadequacy and the displayed SpO₂ or pulse rate value is potentially incorrect.
2. If the screen display "?", it means the signal is unstable, please keep your hands still and retry.
3. PI means Perfusion Index.

Battery Installation


- Slide the battery door cover horizontally along the arrow shown as the picture.
- Install two AAA batteries into the battery compartment. Match the plus (+) and minus (-) signs in the compartment. If the polarities are not matched, damage may cause to the oximeter.
- Close the battery door.

Notes:
1. Please remove the batteries if the pulse oximeter will not be used for long periods of time.
2. Please replace the battery when the power indicator starting flickering.

Using the Lanyard

- Thread thinner end of the lanyard through the hanging hole.
- Thread thicker end of the lanyard through the threaded end before pulling it tightly.

Warnings!
1. Keep the oximeter away from young children. Small items such as the battery door, battery, and lanyard are choking hazards.
2. Do not hang the lanyard from the device's electrical wire.
3. Please notice that the lanyard which is tied to the oximeter may cause strangulation due to excessive length.

Maintenance and Storage

- Replace the batteries in a timely manner when low voltage lamp is lighted.
- Clean surface of the fingertip oximeter before it is used in diagnosis for patients.
- Remove the batteries if the oximeter is not operated for a long time.
- It is best to store the product in -25°C ~ +70°C and ≤ 93% humidity.

5. Keep in dry place. Extreme moisture may affect oximeter lifetime and may cause damage.
6. Disposal of batteries: properly follow any applicable local battery disposal laws.

Cleaning the fingertip pulse oximeter

Please use medical alcohol to clean the silicone touching the finger inside of oximeter with a soft cloth dampened with 70% isopropyl alcohol. Also clean the testing finger using alcohol before and after each test.

Do not pour or spray liquids onto the oximeter, and do not allow any liquid to enter any openings in the device. Allow the oximeter to dry thoroughly before reuse.

The fingertip pulse oximeter requires no routine calibration or maintenance other than replacement of batteries.

The applied products touching the patient's body are required to be disinfected once after each use. The recommended disinfectants include: ethanol 70%, isopropanol 70%, glutaraldehyde-type 2% liquid disinfectants.

Disinfection may cause damage to the equipment and is therefore not recommended for this pulse oximeter unless otherwise indicated in your hospital's servicing schedule. Clean the pulse oximeter before disinfecting it.

CAUTION: Never use ETO or formaldehyde for disinfection.

The use life of the device is five years when it is used for 15 measurements every day and 10 minutes per one measurement. Stop using and contact local

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