
**MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
MANUAL DE USO Y MANTENIMIENTO**

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


DESCRIPTION

ASSEMBLY

1. Insert on base (1) (fig. 2) the stand (2) fixing it with the screw (3) (use the wrench standard supplied).

2. Assemble the optical body (4) to the stand screwing in the small grip handle (5).

3. Connect the plug (8) in the socket (9).

4. Verify that the mains voltage is the same prescribed in the Specifications.

5. Connect the power cable (6) to the power unit (7) and therefore to the mains supply.

USING

1. By (8) knob fig.1, set the wished vision height.
2. Regulate the pan-tilt resistance with the (6-7) knobs, fig 1.
3. Turn "ON" (10) fig.1 the power unit.
4. Using a plane surface object to observe, set the colposcope at an approximate focus position (250 mm). Find the best focus image with focusing knob (5) fig.1.
5. Set the correct interpupillary distance removing or approaching the eyepieces with both the hands until the two images overlap in one circle.
6. Diopter compensation: focus on the object looking through the left eyepiece. Looking through the right eyepiece focusig the object turning the diopter compensation ring (2) fig.1.

MAINTENANCE-CARE

Replacing lamp:

1. Unscrew the protection filter.
2. Pull out the broken lamp.

Attention: Before touch the lamp with the hand, let cool it !

3. When insert the new lamp be carefull to the electrical connection (3) is as shown in figure.

Warning: Use only lamps complying with the original. See Specifications. Never touch the bulb inside with your fingers!

4. Screw the protection filter.


fig.3

Replacing fuses:

Pull out the power cord from the power unit and open the fuses shell as shown in fig.4.

Warning: Use only fuses complying with the original. See Specifications.


fig.4

Care

The quality of the image depends mostly on the cleanliness of the eyepieces and of the front protection filter.


To remove dirt from these, use a soft paintbrush. For persistent spots it is advised to use a wedge of cotton wool soaked in alcohol.

The mechanical parts can be cleaned with up with a soft paint-brush. Keep the colposcope far from corrosive substances.

Always keep the colposcope covered with its dust cover when not in use.

The quality of the materials employed in the construction will guarantee the best performances for several years. If Your instrument does not work correctly, don't hesitate to contact Your GIMA agent.

ELECTRICAL PARTS

SPECIFICATIONS

Eyepieces	10x - 15x
Magnification	6,2 - 9,3
Field of view	35mm - 24mm
Working distance	250mm
Interpupillary distance	da 48mm a 77mm
Diopter compensation	±5 on right eyepieces
Height of vision	da 73 a 104cm
Fine focus	36mm
Filter	Green
Lamp	Halogen 12W, 12V GE M64-FTA 30759
Fuses	T200mA (x2)
Classification	I tipo B
Complayng with	IEC-601-1 93/42 CEE
Main supply	230 V, 50/60 HZ (115 VAC on request)
Weight	5,2kg
Marcatura	CE

Optionals

Eyepieces 15X (pair).

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.


The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.



The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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