

Product Data Sheet



Product name: MULTIPAD Standard
Reference number: 1012-E

Intended use	Universal Patient Positioning System for CT and MR Imaging	
Application site	Head and extremities	
Directives & classification	CE according to MDD – Directive 93/42/ECC Class 1 product	
Applied quality assurance	SN EN ISO 13485:2016	
Instructions for use	Enclosed to the product	
Restrictions of use	Do not use the product with children below the age of 3. Do not use the product for purposes other than those described.	
General safety notes	Non sterile. Do not sterilize. Do not allow direct contact with the skin. Use single-use sanitary cover. Do not use in case of acute open wounds. Do not inflate with compressed air. Position the pump far away from the field of view. Caution: Federal law restricts this device to sale by or on the order of a licensed health-care professional (Rx only).	
Dimensions	Area: 170 mm x 65 mm Variable thickness: 20-50 mm	
Materials	Material TPU EPS Nylon PVC LDPE PP TPE IR LDPE	Component foils, tube pearls Y connection piece pump ball valve body valve lid valve sealing element slit valve packaging
Marking	CE	
Pieces	1	
Packaging	Seal bag	
Packaging dimensions	Length: 470 mm Width: 280 mm Height: 60 mm	
Weight	0.18 kg	
Storage	Store in dry place / room temperature	

Figures are subject to usual tolerances. All information to the best of our current knowledge. They do not guarantee a warranted quality and are no foundation for a contractual relationship. Errors excepted.